

Instructions for use - Servo-Sil 70**Manufacturer**

Servo-Dental GmbH & Co. KG
Rohrstraße 30
D-58093 Hagen (Germany)

Product Description

Servo-Sil 70 is an addition-curing 1:1 kneading silicone for universal use in the dental laboratory.

Indication

For production of silicone indexes and repair models as well as for blocking out undercuts.

Processing

Two equal amounts of component A and B are mixed in order to get a homogeneous mass.

Technical data

Consistency	kneadable
Colour	A: white; B: violet
Mixing ratio	1:1
Mixing time	30 sec.
Processing time	approx. 2 min. (22°C)
Demoulding	10 min (22°C)
Contraction	< 0,16%
Shore hardness	70
Boil-proof	√

Storage

5°C – 25°C

Shelf-life

2 years from date of manufacture, see package.

Available packaging units

REF. 2000600 Kneading silicone Servo-Sil 70, Set: Compound A & B (1 kg each)
REF. 2000605 Kneading silicone Servo-Sil 70, Set: Compound A & B (5 kg each)
REF. 2000610 Kneading silicone Servo-Sil 70, Set: Compound A & B (7.5 kg each)

25.03.2014

Rev.1