SD-Snap Riegel

SD-Snap - in latch

SD-Snap Riegel

SD-Snap - in latch

Indikation:

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Indication:

The SD-SNAP-IN LATCH is an extracoronaral retention, supporting and splinting unit for use in partial dentures and implant-borne superstructures.

Particular advantages:


The SD-SNAP-IN extracoronaral latch locks into place automatically. Light pressure on the pin is sufficient to unlock it automatically and it can be removed without having to touch the latch.

The primary component is made of fully combustible plastic, high fusing for casting to precious and semi-precious alloys or platinum/iridium alloy for all metals including titanium. The secondary component is made of titanium and is designed for adhering into a chrome cobalt denture base.

Two versions of the SD-SNAP-IN LATCH are available, with a right-hand pin (red) and a left-hand pin (blue).

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
Tooth guide to help you choose and order the right colour of the SD-Snap in latch, depending on the position of the axle.


Fig. 1

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Waxing-up:

The crown is waxed-up using standard procedures and is milled circumferentially.

Fig. 2

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The angle of insertion of the SD-SNAP-IN LATCH is first established using the paralleling mandrel (Order No. 8453).

Fig. 3

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The primary component is waxed onto the model. The backplate must be blocked-out properly with wax. To ensure that the periodontal tissues are not impinged on, the backplate can be extended at an angle of 45° (Order No. 8451/0 red or 8452/0 blue).

Fig. 4

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Investing:

The crowns with the primary component are invested and cast using standard methods. A ceramic pin (Order No. 8458) is inserted into the aperture in the primary component to ensure that the pin aperture is cast properly and that precise axial guidance is obtained.


Fig. 5

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Casting:

Devest the crowns after casting. Carefully remove any casting defects and sandblast the primary component with fine grain abrasive or glass beads. Use a reamer (Order No. 8454).

Fig. 6

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Trimming and Finishing:

After applying the facings and finishing the primary, replace the crowns on the model and prepare it for duplication. Wax is applied beneath the primary component only.

Fig.7

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Preparing for Duplication:

The brass duplicating aid (Order No. 8452/4) is placed on the secondary component and the reamer (Order No.8454) is used to check that the secondary component corresponds to the external duplicating component. If necessary, it should be adjusted. The ledge on the bottom of the duplicating aid acts as a spacer for the latch pin safety wire. That side of the duplicating aid facing the alveolar ridge should be blocked out with wax.

Fig. 8

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Silicone is well suited as a duplicating material. The duplicating pin should only protrude from the duplicating aid by about 1mm.

Fig. 9

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Refractory Model:

The area where the latch pin is pressed out can be seen clearly on the investment model.

Fig. 10

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Wax Pattern:

The pattern is now waxed-up on the investment model, ensuring that the wax does not obstruct the pin aperture. Please note that shear distribution arms and milled interlocks should be used for free-end cases.


Fig. 11

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The denture base pattern is then invested and cast according to the manufacturer's recommendations. The duplicating aid is first removed and all the wax carefully cleaned off the crowns with the primary component. The crowns with the primary component are placed in position and the actual secondary component placed on it.

Fig. 12

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Adhering:

In order to prepare for adhering with SD-MASTERBOND (Order No. 8300), that area of the denture base to receive the secondary component should be sandblasted (250 µm). The plastic sleeve (Order No. 8457) should now be pressed onto the latch pin. The basal end of the retention pin (Order No. 8459) should be bent. The clean receptacle for the denture base is coated with SD-MASTERBOND and the denture base placed on the secondary component.

Fig. 13

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Finishing:

Remove any excess SD-MASTER-BOND. Finish the denture using standard methods. The plastic sleeve around the latch pin ensures that it moves smoothly (if necessary, trim it slightly with a reamer). The pin can be shortened as required.

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Starter set:

Order No.: 8450
 Contents: 2 SD-SNAP-IN LATCHES, red
 2 SD-SNAP-IN LATCHES, blue
 1 paralleling mandrel
 1 reamer
 1 duplicating aid, red
 1 duplicating aid, blue
 4 ceramic pins

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Separate components:

SD-SNAP-IN LATCH, red	Order No.: 8451
SD-SNAP-IN LATCH, blue	Order No.: 8452
Primary part, red	Order No.: 8451/3
Primary part, blue	Order No.: 8452/3
Primary part, HFA - red	Order No.: 8451/5
Primary part, HFA - blue	Order No.: 8452/5
Primary part, Pt./Ir. - red	Order No.: 8451/6
Primary part, Pt./Ir. - blue	Order No.: 8452/6

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
Best-Nr.: 8451 /4
 Best-Nr.: 8452 /4
 Best-Nr.: 8460

Duplicating aid, red
Duplicating aid, blue
Duplicating pin
(for red and blue)

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Secondary part, red
Primary part, red
Primary part extension, red
Leaf spring
Spiral spring
Pin, red
Plastic sleeve

Order No.:8451/1
 Order No.: 8451/3
 Order No.: 8451/0
 Order No.: 8455
 Order No.: 8456
 Order No.: 8451/2
 Order No.: 8457

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Secondary part, blue
Primary part, blue
Primary part extension, blue
Leaf spring
Spiral spring
Pin, blue
Plastic sleeve

Order No.: 8452/1
 Order No.: 8452/3
 Order No.:8452/0
 Order No.: 8455
 Order No.: 8456
 Order No.: 8452/2
 Order No.: 8457

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Please note:

If necessary, both the pin and spring in the secondary component can be replaced. In order to do so, remove the safety wire from the underside so that the latch pin can be withdrawn. To remove the spring, first lever it upwards from the interdental side with a No. 531 M instrument and press the spring downwards and out. To insert the spring, introduce it from beneath, until it locks in place automatically.

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Servo-Dental GmbH & Co. KG