

SD-TeleC**lick**

(nach Ztm. Peter Hackstein)
(Idea from MDT Peter Hackstein)

Maßstab 1:1

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Natural scale

-  Spring: 2.4 mm high, 2.0 mm wide
-  Female: 2.7 mm high, 2.4 mm wide, 0.5mm deep
-  Male: 2.7 mm high, 2.4 mm wide, 0.5mm deep

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Indication

The SD-TeleClick is an intracoronal retaining element for telescopic crowns. It is also suitable on the side of a bar in combination with implants. The SD-TeleClick for precious and semi precious alloys is made of HFA and for non precious alloys it is made of Pt/Ir.

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Advantages:

- Extremely small dimensions of the assembly
- Confident patient due to the click that secures the denture in the final position
- Exchanging the spring is made easy by using the spring remover
Can be used with nearly every alloy on the market
- For high and semi precious alloys using HFA male and female
- For non precious using platinum-iridium male and female
- The spring is made of 22 carat gold and is only 0.2mm thick

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
Fig. 1 & 2:

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The SD-TeleClick female is positioned using the parallelling tool and waxed into the telescopic / conus crown or into a bar.

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The upper edge of the female must remain free of wax (see fig. 1).


Fig. 3:

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Crowns or sub-structure are sprued and invested as usual.

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Prior to investing wetting agent must be applied to the female and dried.

Fig. 4:

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The cast object.

Please note: The female is only 0.5mm thin.

The female should only be sandblasted very carefully with 50µ glass beads / max. 2 bar (30 psi).

The occlusive ridge and the female itself is not allowed to be harmed due to milling or polishing.

Fill the female with a little drop of resin and after finishing remove by heating slightly.


Fig. 5 & 6:

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The male is provided with a try-in spring and is positioned into the female. The try-in spring guarantees the precise positioning of the male and the modulation can be removed easily.

Important: *The occlusive edge of the male has to be waxed up by 0.5mm providing enough space to use the spring remover.*

The SD-TeleClick is also available as a resilient version. The female is 3.3mm high, the male is fitted with a red try-in spring and is also 3.3mm high.

Fig. 7 & 8:

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After lifting off the secondary structure remove the try-in spring using the spring remover. Be sure that there is no wax left in the area for spring removal. Sprue, invest and cast as usual. For the wax try-in use the black try-in spring made of plastics instead of the original gold spring because this is very difficult to remove without distorting the wax set-up.



Fig. 9 & 10

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After finishing the job the gold spring is pushed over the male using the positioning aid provided with the spring.

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If the spring has to be exchanged use the spring remover as shown.



Fig. 15:

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Insert the new spring with the positioning aid.

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- | | |
|--|------------------------|
| 1. SD-TeleClick Matrize resi., Pt./Ir. | Best.-Nr.: 138 121 200 |
| 2. SD-TeleClick Matrize resi., HSL | Best.-Nr.: 138 121 300 |
| 3. SD-TeleClick Matrize starr, Pt./Ir. | Best.-Nr.: 138 131 200 |
| 4. SD-TeleClick Matrize starr, HSL | Best.-Nr.: 138 131 300 |
| 5. SD-TeleClick Patrize Pt./Ir. | Best.-Nr.: 138 141 200 |
| 6. SD-TeleClick Patrize HSL | Best.-Nr.: 138 141 300 |
| 7. SD-TeleClick Goldfeder | Best.-Nr.: 138 784 200 |
| 8. SD-TeleClick Parallelhalter | Best.-Nr.: 138 593 100 |
| 9. SD-TeleClick Federentferner | Best.-Nr.: 138 953 100 |

- | | |
|---------------------------------------|------------------------|
| 1. SD-TeleClick Female resi., Pt./Ir. | Order No.: 138 121 200 |
| 2. SD-TeleClick Female resi., HFA | Order No.: 138 121 300 |
| 3. SD-TeleClick Female rigid, Pt./Ir. | Order No.: 138 131 200 |
| 4. SD-TeleClick Female rigid, HFA | Order No.: 138 131 300 |
| 5. SD-TeleClick Male Pt./Ir. | Order No.: 138 141 200 |
| 6. SD-TeleClick Male HFA | Order No.: 138 141 300 |
| 7. SD-TeleClick Gold spring | Order No.: 138 784 200 |
| 8. SD-TeleClick Paralleling tool | Order No.: 138 593 100 |
| 9. SD-TeleClick Spring remover | Order No.: 138 953 100 |

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SD-TeleClick standard

Order No.: 138 011 200

Starter set SD-TeleClick, rigid, Pt./Ir.,
castable to NP and Titanium,
consisting of 4x SD-TeleClick complete,
one paralleling tool and one spring remover.

Order No.: 138 021 200

Starter set SD-TeleClick, resi, Pt./Ir.,
castable to NP and Titanium,
consisting of 4x SD-TeleClick complete,
one paralleling tool and one spring remover.

Order No.: 138 031 300

Starter set SD-TeleClick, rigid, HFA,
castable to precious and semi precious,
consisting of 4x SD-TeleClick complete,
one paralleling tool and one spring remover.

Order No.: 138 041 300

Starter set SD-TeleClick, resi, HFA,
castable to precious and semi precious,
consisting of 4x SD-TeleClick complete,
one paralleling tool and one spring remover.


	SD-TeleClick rigid, Pt./Ir. single item	Order No.: 138 501 200
	SD-TeleClick rigid, HFA single item	Order No.: 138 501 300
	SD-TeleClick resi., Pt./Ir. single item	Order No.: 138 511 200
	SD-TeleClick resi., HFA single item	Order No.: 138 511 300

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Servo-Dental GmbH & Co. KG