








































Ball-Anchors horizontal: normal Ø 2,5mm, micro Ø 1,8mm

Item No.		Description
2.58.15.98.03		<p><b>Introductory Kit Laboratory</b></p> <ul style="list-style-type: none"> <li>16 Castable Ball-Anchor Pivots , all sizes (normal/mikro)</li> <li>2 Castable bars regular size (1 straight/1 bevelled)</li> <li>2 Castable bars micro size (1 straight/1 bevelled)</li> <li>2 Basic putting profiles normal (1 upper/1 lower)</li> <li>2 Basic putting profiles micro (1 upper/1 lower)</li> <li>1 Basic profile "Monoblock" normal</li> <li>1 Basic profile "Monoblock" micro</li> <li>2 Matrix cages "single piece", plastic, normal</li> <li>2 Matrix cages "single piece", plastic, mikro</li> <li>2 Ball Anchors-horizontal "single", normal</li> <li>2 Ball Anchors-horizontal "single", mikro</li> <li>4 plastic positioning rings (2 regular size/2 mikro)</li> <li>4 BA-horizontal white caps (standard retention) (2 normal/2 mikro)</li> <li>4 BA-horizontal pink caps (soft retention) (2 normal/2 mikro)</li> <li>4 BA-horizontal yellow caps (medium retention) (2 normal/2 mikro)</li> <li>4 BA-horizontal green caps (resilient) (2 normal/2 mikro)</li> <li>4 BA-horizontal working caps (black)(2 normal/2 mikro)</li> <li>2 Stainless steel housings normal (for soldering and for resin)</li> <li>2 Stainless steel housings mikro (for soldering and for resin)</li> <li>4 Aluminium spacers for BA-horizontal bar (2 normal/ 2 mikro)</li> <li>4 Ball-Anchors vertical (2 normal/ 2 long)</li> <li>2 Elongation for Ball-Anchors vertical</li> <li>4 BA-vertical white caps for duplication (standard retention)</li> <li>2 analogue pivots for stone cast for BA-horizontal normal</li> <li>2 analogue pivots for stone cast for BA-horizontal mikro</li> <li>2 analogue pivots for stone cast for BA-vertical</li> <li>3 Keys for parallelometer (horizontal normal &amp; mikro / vertical)</li> <li>3 Reamers for regulating caps (horizontal normal &amp; mikro / vertikal)</li> <li>2 Tools for inserting caps (horizontal normal &amp; mikro / vertical)</li> <li>4 Aluminium disks for Ball-Anchor "single-type"</li> <li>1 Tweezers</li> </ul>
2.41.03.98.03		<p><b>BA-Kit horizontal, normal, single parts</b></p> <ul style="list-style-type: none"> <li>4 BA-horizontal singles normal size</li> <li>4 BA-caps horizontal normal, pink - soft</li> <li>4 Stainless steel housings normal (for soldering and for resin)</li> <li>4 Plastic positioning rings normal</li> </ul>
2.41.04.98.03		<p><b>BA-Kit horizontal, micro, single parts</b></p> <ul style="list-style-type: none"> <li>4 BA-horizontal singles micro size</li> <li>4 BA-caps horizontal micro, pink - soft</li> <li>4 Stainless steel housings mikro (for soldering and for resin)</li> <li>4 Plastic positioning rings mikro</li> </ul>
2.41.05.98.03		<p><b>BA-Eco-Kit horizontal, normal</b></p> <ul style="list-style-type: none"> <li>1 Castable bar normal, straight</li> <li>1 Castable bar normal, bevelled</li> <li>4 BA-horizontal,normal, white caps (standard retention)</li> </ul>

2.41.06.98.03		<b>BA-Eco-Kit horizontal, micro</b> 1 Castable bar micro, straight 1 Castable bar micro, bevelled 4 BA-horizontal, micro, white caps (standard retention)
2.41.07.98.03		<b>BA-Kit horizontal, with matrix cage "single piece", normal</b> 1 Castable bar normal, straight 1 Castable bar normal, bevelled 4 Matrix cages "single piece", plastic, normal 4 Plastic positioning rings normal 4 BA-horizontal,normal, white caps (standard retention)
2.41.08.98.03		<b>BA-Kit horizontal, with matrix cage "single piece", micro</b> 1 Castable bar micro, straight 1 Castable bar micro, bevelled 4 Matrix cages "single piece", plastic, micro 4 Plastic positioning rings normal 4 BA-horizontal,micro, white caps (standard retention)
2.41.15.98.03		<b>BA-horizontal-Kit for caps horizontal, normal, assorted</b> 6 pieces BA-horizontal caps normal per colour: white-standard, pink-soft,yellow-medium, green-resilient
2.41.16.98.03		<b>BA-horizontal-Kit for caps horizontal, micro, assorted</b> 6 pieces BA-horizontal caps micro per colour: white-standard, pink-soft,yellow-medium, green-resilient
2.40.01.20.03		2 BA-horizontal stainless steel housings, for soldering and for resin, normal
2.40.02.20.03		2 BA-horizontal stainless steel housings, for soldering and for resin, micro
2.40.19.75.03		2 BA horizontal, Castable bar normal, bevelled, plastic

2.40.20.75.03		2 BA horizontal, Castable bar normal, straight, plastic
2.40.22.75.03		2 BA horizontal, Castable bar micro, bevelled, plastic
2.40.23.75.03		2 BA horizontal, Castable bar micro, straight, plastic
2.40.21.75.03		4 BA-horizontal singles normal size, plastic, green
2.40.89.75.03		4 BA-horizontal singles micro size, plastic, red
2.40.26.75.03		6 BA-horizontal,normal, pink caps (soft retention)
2.40.27.75.03		6 BA-horizontal,normal, yellow caps (medium retention)
2.40.28.75.03		6 BA-horizontal,normal, white caps (standard retention)
2.40.29.75.03		6 BA-horizontal,normal, green caps (resilient)
2.40.30.75.03		6 BA-horizontal,micro, pink caps (soft retention)
2.40.31.75.03		6 BA-horizontal,micro, yellow caps (medium retention)
2.40.32.75.03		6 BA-horizontal,micro, white caps (standard retention)
2.40.33.75.03		6 BA-horizontal,micro, green caps (resilient)
2.40.47.38.03		10 Aluminium spacers for BA-horizontal bar normal
2.40.48.38.03		10 Aluminium spacers for BA-horizontal bar micro
2.40.49.38.03		30 Disks for Ball-Anchor "single-type"
2.40.50.20.03		2 analogue pivots for stone cast for BA-horizontal normal
2.40.51.20.03		2 analogue pivots for stone cast for BA-horizontal micro
2.40.53.75.03		6 BA-horizontal working caps (black) normal -only for laboratory use-
2.40.54.75.03		6 BA-horizontal working caps (black) micro -only for laboratory use-
2.40.57.32.03		Reamer for regulating caps, normal
2.40.58.32.03		Reamer for regulating caps, micro
2.40.60.98.03		Tool for inserting caps, horizontal, normal und micro
2.40.62.98.03		Key for parallelometer, horizontal, normal
2.40.63.98.03		Key for parallelometer, horizontal, micro
2.40.65.32.03		Tool for testing caps retention, horizontal, normal
2.40.66.32.03		Tool for testing caps retention, horizontal, micro
2.41.25.98.03	<b>TOOLKIT</b>	<b>Tool-Kit for the dentist (5 tools)</b> Contains: 2.40.57.32.03, 2.40.58.32.03, 2.40.60.98.03, 2.40.61.98.03, 2.40.67.75.03